



Recombinant Human Beta-APP39 Protein [His, GST] (DAG-WT3491)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	His, GST
Applications	N/A
Molecular Weight	34 kDa
Format	Liquid
Size	100 µg
Buffer	50 mM Tris, 100 mM NaCl, 20% glycerol, 0.05% Tween 20, pH 9.5
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Amyloid-beta precursor protein; APP; APP-beta; Beta-APP39
----------	-----------------------------------------------------------

GENE INFORMATION

UniProt ID	P05067-1
------------	----------