



Recombinant Human Beta-APP40 Protein [His, GST] (DAG-WT3492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 80% as determined by SDS-PAGE
Conjugate	His, GST
Applications	N/A
Molecular Weight	33 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile 50mM Tris, 500mM NaCl, pH 7.5
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Amyloid-beta precursor protein; APP; APP-beta; Beta-APP40
----------	---

GENE INFORMATION

UniProt ID	P05067-1
------------	----------