



Recombinant Human Beta-APP42 Protein [His, GST] (DAG-WT3493)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 80% as determined by SDS-PAGE
Conjugate	His, GST
Applications	N/A
Molecular Weight	34 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile PBS, 10% glycerol, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Amyloid-beta precursor protein; APP; APP-beta; Beta-APP42
----------	---

GENE INFORMATION

UniProt ID	P05067-1
------------	----------