



Recombinant Human Beta-APP Protein [hFc] (DAG-WT3494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	hFc
Applications	N/A
Molecular Weight	155 kDa
Format	Lyophilized
Size	50 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Amyloid-beta precursor protein; APP; APP-beta
----------	---

GENE INFORMATION

UniProt ID	NP_958816.1
------------	-------------