



Recombinant Human Clusterin Protein (CLU34) [His] (DAG-WT3495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	38 kDa
Format	Lyophilized
Size	50 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	CLU; Clusterin
----------	----------------

GENE INFORMATION

UniProt ID	NP_001822.2
------------	-------------