



Recombinant Human BACE1 Protein [hFc] (DAG-WT3496)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	hFc
Applications	N/A
Molecular Weight	101 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	BACE1; Beta-secretase 1
----------	-------------------------

GENE INFORMATION

UniProt ID	NP_036236.1
------------	-------------