



# Recombinant Human BACE1 Protein (DAG-WT3498)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	65 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 8.0
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	BACE1; Beta-secretase 1
----------	-------------------------

## GENE INFORMATION

UniProt ID	NP_036236.1
------------	-------------