



Biotinylated Recombinant Human ApoE3 Protein [Avi, His] (DAG-WT3500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 75% as determined by SDS-PAGE
Conjugate	Avi, His
Applications	N/A
Molecular Weight	34 kDa
Format	Lyophilized
Size	20 µg, 100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Apolipoprotein E; APOE
----------	------------------------

GENE INFORMATION

UniProt ID	P02649
------------	--------