



Recombinant Mouse BDNF Protein [His] (DAG-WT3501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Mouse
Purity	> 75% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	51.3 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Brain derived neurotrophic factor; BDNF
----------	---

GENE INFORMATION

UniProt ID	EDL27780.1
------------	------------