



Recombinant Mouse RAGE Protein [His] (DAG-WT3502)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Mouse
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	27.4 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Advanced glycosylation end product-specific receptor isoform a precursor; RAGE
----------	--

GENE INFORMATION

UniProt ID	NP_031451.2
------------	-------------