



# Recombinant Human RAGE Protein [His] (DAG-WT3503)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	50 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	Advanced glycosylation end product-specific receptor isoform a precursor; RAGE
----------	--

## GENE INFORMATION

UniProt ID	NP_001127
------------	-----------