



Recombinant Human RAGE Protein [hFc] (DAG-WT3504)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 70% as determined by SDS-PAGE
Conjugate	hFc
Applications	N/A
Molecular Weight	85 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Advanced glycosylation end product-specific receptor isoform a precursor; RAGE

GENE INFORMATION

UniProt ID NP_001127