



Biotinylated Recombinant Human RAGE Protein [Avi, hFc] (DAG-WT3506)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	Avi, hFc
Applications	N/A
Molecular Weight	85 kDa
Format	Lyophilized
Size	20 µg, 100 µg
Buffer	Lyophilized from sterile PBS, pH 7.
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Advanced glycosylation end product-specific receptor isoform a precursor; RAGE
----------	--

GENE INFORMATION

UniProt ID	NP_001127
------------	-----------