



Recombinant Mouse AChE Protein [His] (DAG-WT3507)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Mouse
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	66.2 kDa
Format	Lyophilized
Size	20 µg
Buffer	Lyophilized from sterile PBS, pH 7.
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Acetylcholinesterase precursor; AChE
----------	--------------------------------------

GENE INFORMATION

UniProt ID	NP_033729.1
------------	-------------