



Recombinant Human SOD1 (DAG-WT3508)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	17 kDa
Bio-activity	≥10,380 IU/mg
Format	Lyophilized
Size	1 mg, 10 mg, 100 mg
Buffer	Lyophilized from sterile water
Storage	RT

BACKGROUND

Keywords	SOD1; Superoxide dismutase 1
----------	------------------------------