



Recombinant Alkali-tolerant Protein A (DAG-WT3509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	28 kDa
Format	Lyophilized
Size	1 mg, 10 mg, 100 mg
Buffer	Lyophilized from 5 mM PB
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Protein A; SpA; Alkali-tolerant protein A
----------	---