



# Recombinant Core Streptavidin (DAG-WT3510)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	13.5 kDa
Bio-activity	≥15 U/mg
Format	Lyophilized
Size	1 mg, 10 mg, 100 mg
Buffer	Lyophilized from 5 mM PB
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	Streptavidin; SA
----------	------------------