



Recombinant Bovine Enterokinase [His] (DAG-WT3512)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	30 kDa
Bio-activity	≥40000 U/mg
Format	Lyophilized
Size	200 U, 1000 U
Buffer	Lyophilized from 50 mM Tris-HCl, pH 8.0
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Bovine enterokinase; Enterokinase; ENTK
----------	---