



Recombinant Human Fibroblast Growth Factor (FGF1) (DAG-WT3514)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Molecular Weight	16 kDa
Bio-activity	EC50 ≤0.5 ng/ml
Format	Lyophilized
Size	10 µg, 100 µg, 1 mg
Buffer	Lyophilized from 10 mM PBS
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Human fibroblast growth factor; FGF1
----------	--------------------------------------