



Recombinant *Anaplasma phagocytophilum* Chaperonin GroEL [His, Myc] (DAG-WT3518)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His, Myc
Applications	ELISA
Molecular Weight	28.6 kDa
Format	Lyophilized
Size	20 µg, 100 µg, 1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords *Anaplasma phagocytophilum*; Chaperonin GroEL; groEL; 60 kDa chaperonin; Chaperonin-60; Cpn60

GENE INFORMATION

UniProt ID Q2GL94