



Recombinant *Anaplasma phagocytophilum* Large ribosomal subunit protein uL22 [His, Myc] (DAG-WT3519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity > 90% as determined by SDS-PAGE

Conjugate His, Myc

Applications ELISA

Molecular Weight 17.2 kDa

Format Lyophilized

Size 20 µg, 100 µg, 1 mg

Buffer Tris-based buffer, 50% glycerol

Storage Store at -20°C to -80°C

BACKGROUND

Keywords *Anaplasma phagocytophilum*; Large ribosomal subunit protein uL22; rplV; 50S ribosomal protein L22

GENE INFORMATION

UniProt ID Q2GL54