



Recombinant Hepatitis B virus Surface Antigen L-protein (genotype D) (DAG-WT3522)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, WB
Reconstitution	Added 500 µL of water to the vial that makes a antigen solution at 0.2 mg/mL in PBS
Format	Lyophilized
Size	100 µg
Buffer	PBS, pH 7.3, 1% sucrose
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	HBsAg L-protein; HBsAg; Large HBV surface protein
----------	---