



Recombinant Hepatitis B virus Surface Antigen S-protein (adw) (DAG-WT3523)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	10mM sodium phosphate, 140mM NaCl, pH 6.4
Preservative	0.1% sodium azide
Storage	Store at 4°C

BACKGROUND

Keywords	HBsAg S-protein; HBsAg; Small HBV surface protein
----------	---