



Recombinant Human Smooth Muscle Actin Protein (ACTA2) [His] (DAG-WT3525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------|---|
| Purity | > 95% as determined by SDS-PAGE |
| Conjugate | His |
| Applications | N/A |
| Molecular Weight | 45.5 kDa |
| Format | Lyophilized powder |
| Size | 50 µg, 200 µg, 1 mg |
| Buffer | 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose |
| Preservative | Proclin300 |
| Storage | Store at -20°C |

BACKGROUND

Keywords Smooth muscle aortic alpha-actin; ACTA2

GENE INFORMATION

UniProt ID P62736