



Testosterone BSA Conjugate (DAG-WT3527)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	Immunogen; SDS-PAGE; WB
Reconstitution	Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex
Format	Lyophilized
Size	50 µg, 200 µg, 1 mg
Buffer	PBS, pH7.4
Preservative	None
Storage	Short-term (2 weeks) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles

BACKGROUND

Keywords	Testosterone; 17β-Hydroxyandrost-4-en-3-one; Androst-4-en-17β-ol-3-one
----------	--