



Recombinant Human UCHL1 Protein (DAG-WT3528)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------|---------------------------------|
| Purity | > 85% as determined by SDS-PAGE |
| Conjugate | TBD |
| Applications | ELISA |
| Format | Liquid |
| Size | 100 µg |
| Buffer | 10 mM PBS, pH 7.4 |
| Preservative | None |
| Storage | Store at -20°C |

BACKGROUND

| | |
|----------|--|
| Keywords | UCHL1; Ubiquitin carboxyl-terminal hydrolase isozyme L1; PGP 9.5 |
|----------|--|