



Recombinant HIV gp41+gp36 antigen (Detection) (DAG-WT3529)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	TBD
Applications	LFIA (Detection)
Format	Liquid
Size	1 mg
Buffer	10 mM PBS, pH 7.4
Preservative	None
Storage	Short-term (2 weeks) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles

BACKGROUND

Keywords	HIV; Envelope glycoprotein; gp41; gp36
----------	--