



# Recombinant Glycine max Gly m 3 (DAG-WT3539)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Conjugate</b>	TBD
<b>Applications</b>	ELISA, CLIA, LFIA WB
<b>Molecular Weight</b>	15.2 kDa
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS, pH 7.3
<b>Preservative</b>	0.1% ProClin300
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

<b>Keywords</b>	Gly m 3; Glycine max; Profilin
-----------------	--------------------------------