



UR-144 [HRP] (DAG-WT354)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	HRP
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	UR-144 is a drug invented by Abbott Laboratories, that acts as a selective full agonist of the peripheral cannabinoid receptor CB2, but with much lower affinity for the psychoactive CB1 receptor.
Keywords	UR-144; TMCP-018; KM-X1; MN-001; YX-17