



# Recombinant Phaseolus vulgaris Pha v 3 (DAG-WT3541)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	20.7 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.3
Preservative	0.1% ProClin300
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	Pha v 3; Phaseolus vulgaris; Non-specific lipid transfer protein type 1
----------	---