



# Recombinant *Pyrus communis* Pyr c 1 (DAG-WT3543)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90% as determined by SDS-PAGE
---------------	---------------------------------

<b>Conjugate</b>	TBD
------------------	-----

<b>Applications</b>	ELISA, CLIA, LFIA WB
---------------------	----------------------

<b>Molecular Weight</b>	18.6 kDa
-------------------------	----------

<b>Format</b>	Liquid
---------------	--------

<b>Size</b>	1 mg
-------------	------

<b>Buffer</b>	PBS, pH 7.3
---------------	-------------

<b>Preservative</b>	0.1% ProClin300
---------------------	-----------------

<b>Storage</b>	Store at -20°C to -80°C
----------------	-------------------------

## BACKGROUND

<b>Keywords</b>	Pyr c 1; <i>Pyrus communis</i> ; Pathogenesis-related protein
-----------------	---