



Recombinant Cladosporium herbarum Cla h 8 (DAG-WT3547)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	29.5 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.3
Preservative	0.1% ProClin300
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Cla h 8; Cladosporium herbarum; Mannitol dehydrogenase
----------	--