



Recombinant Cynodon dactylon Cyn d 12 (DAG-WT3551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	15.3 kDa
Format	Liquid
Size	1 mg
Buffer	Phosphate buffer saline, pH 8.0, and 8 M Urea
Preservative	0.1% ProClin300
Storage	Store at 2°C to 8°C

BACKGROUND

Keywords	Cyn d 12; Cynodon dactylon; Profilin
----------	--------------------------------------