



Recombinant Human Inhibin β chain B (DAG-WT3553)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------|---------------------------------|
| Purity | > 90% as determined by SDS-PAGE |
| Conjugate | TBD |
| Applications | ELISA, CLIA, LFIA WB |
| Molecular Weight | 16.9 kDa |
| Format | Liquid |
| Size | 1 mg |
| Buffer | PBS, pH 7.3 |
| Preservative | 0.1% ProClin300 |
| Storage | Store at -20°C to -80°C |

BACKGROUND

| | |
|----------|------------------------|
| Keywords | Inhibin; beta B; INHBB |
|----------|------------------------|