



Recombinant Human Inhibin B ($\alpha + \beta$ chain B) (DAG-WT3554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	34.2 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.3
Preservative	0.1% ProClin300
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Inhibin; beta B; INHBB
----------	------------------------