



Recombinant Deamidated gliadin peptide (20 kDa) (DAG-WT3559)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	19.9 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.3
Preservative	0.1% ProClin300
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Deamidated gliadin peptide; DGP; Celiac disease
----------	---