



Recombinant Deamidated gliadin peptide (33 kDa) (DAG-WT3560)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA, CLIA, LFIA WB
Molecular Weight	33.2 kDa
Format	Liquid
Size	1 mg
Buffer	Tris-buffered saline, pH 8.0
Preservative	0.1% ProClin300
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Deamidated gliadin peptide; DGP; Celiac disease
----------	-------------------------------------------------