



Recombinant Heartland virus Nucleoprotein [His] (DAG-WT3562)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA, Immunogen
Molecular Weight	25 kDa
Format	Liquid
Size	100 µg, 1 mg
Buffer	50mM Tris-HCl, pH8.0, 150mM NaCl and 8M Urea
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Heartland virus; Nucleoprotein; Nucleocapsid protein
----------	--

GENE INFORMATION

UniProt ID	BCU92821
------------	----------