



Recombinant Heartland virus Nucleoprotein [His] (DAG-WT3563)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------|---------------------------------|
| Purity | > 95% as determined by SDS-PAGE |
| Conjugate | His |
| Applications | WB, ELISA, Immunogen |
| Molecular Weight | 27.5 kDa |
| Format | Liquid |
| Size | 100 µg, 1 mg |
| Buffer | PBS |
| Preservative | None |
| Storage | Store at -20°C to -80°C |

BACKGROUND

| | |
|----------|--|
| Keywords | Heartland virus; Nucleoprotein; Nucleocapsid protein |
|----------|--|

GENE INFORMATION

| | |
|------------|----------|
| UniProt ID | BCU92821 |
|------------|----------|