



Recombinant HIV type 1 p31 Integrase [GST] (DAG-WT3565)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	GST
Applications	ELISA, WB
Molecular Weight	57 kDa
Format	Liquid
Size	1 mg
Buffer	50mM Tris-HCl, pH8.0, 150mM NaCl and 8M Urea
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Human immunodeficiency virus; HIV; HIV-1; Integrase; p31
-----------------	--