



Recombinant Heartland virus Glycoprotein G1 [His] (DAG-WT3568)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Format	Liquid
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Heartland virus; Glycoprotein G1; Glycoprotein N
----------	--

GENE INFORMATION

UniProt ID	J3WAX0
------------	--------