



JWH-018 N-pentanoic acid Standard Solution (DAG-WT358)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	ELISA
Molecular Weight	371.4 g/mol
Format	Liquid
Concentration	1 µg/mL
Size	0.5 ml
Buffer	20mM sodium phosphate, 150mM NaCl, pH7.4
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Introduction	JWH 018 is a mildly selective agonist of the peripheral cannabinoid (CB2) receptor, derived from the aminoalkylindole WIN 55,212-2.
Keywords	JWH 018; 1-pentyl-3-(1-naphthoyl)indole; AM678