



Inactivated Coronavirus OC43 Culture Fluid (DAG-WT3616)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ TCID ₅₀ /mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Coronavirus OC43; Human coronavirus; HCoV
----------	---