



# Inactivated Adenovirus type 3 Culture Fluid (DAG-WT3623)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                                    |
|--------------|------------------------------------|
| Conjugate    | N/A                                |
| Applications | Control                            |
| Bio-activity | $\geq 10^6$ TCID <sub>50</sub> /mL |
| Format       | Liquid                             |
| Size         | 1 mL                               |
| Buffer       | DMEM                               |
| Preservative | None                               |
| Storage      | Store at -20°C                     |

## BACKGROUND

|          |                 |
|----------|-----------------|
| Keywords | Adenovirus; ADV |
|----------|-----------------|