



Inactivated *Neisseria gonorrhoeae* Culture Fluid (DAG-WT3640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Neisseria gonorrhoeae</i> ; <i>N. gonorrhoeae</i>
----------	--