



# Inactivated *Bacillus subtilis* Culture Fluid (DAG-WT3677)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Keywords	<i>Bacillus subtilis</i> ; B. subtilis
----------	----------------------------------------