



Inactivated Streptococcus dysgalactiae Culture Fluid (O90R) (DAG-WT3687)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Streptococcus dysgalactiae; S. dysgalactiae
----------	---