



# Inactivated Streptococcus dysgalactiae Culture Fluid (DAG-WT3688)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                    |
|--------------|--------------------|
| Conjugate    | N/A                |
| Applications | Control            |
| Bio-activity | $\geq 10^6$ CFU/mL |
| Format       | Liquid             |
| Size         | 1 mL               |
| Buffer       | DMEM               |
| Preservative | None               |
| Storage      | Store at -20°C     |

## BACKGROUND

|          |                                             |
|----------|---------------------------------------------|
| Keywords | Streptococcus dysgalactiae; S. dysgalactiae |
|----------|---------------------------------------------|