



Inactivated *Corynebacterium glutamicum* Culture Fluid (DAG-WT3691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Corynebacterium glutamicum</i> ; <i>C. glutamicum</i> GS
----------	---