



Inactivated *Ligilactobacillus saerimneri* Culture Fluid (DAG-WT3692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Ligilactobacillus saerimneri</i> ; <i>L. saerimneri</i>
-----------------	--